Handheld Base Unit
Operator's Manual

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Revision History

Manual version	Revision date	Changes
Version 1.0	2008-07	Initial Version
Version 2.0	2015-05	Updated for launch of CoaguChek Pro II for the following languages only: DA, EN, DE, FR, ES, FI, IT, NL, NO, PT, SV; added reference to CoaguChek XS Pro and CoaguChek Pro II systems, updated safety-related information, updated section "Contact Roche", minor changes.
Version 3.0	2017-06	Update: added pictures of new cobas h 232 meter (HBM 4.x platform), updated safety-related information, minor editorial changes.
Version 4.0	2018-09	Update: new hardware, updated cleaning/disinfection section, minor editorial changes
Version 4.1	2021-11	Implementation of IVDR related requirements (IVDR: in vitro diagnostics regulation).

Handheld Base Unit

Operator's Manual

Version 4.1

© 2008-2021 Roche Diagnostics GmbH. All rights reserved.

The contents of this document, including all graphics, are the property of Roche Diagnostics. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Roche Diagnostics.

Roche Diagnostics has made every reasonable effort to ensure that all the information contained in this manual is correct at the time of printing. However, Roche Diagnostics reserves the right to make any changes necessary without notice as part of ongoing product development.

Electronic user documentation can be downloaded using the eLabDoc e-service on Dialog:

www.dialog.roche.com/

For more information, contact your local affiliate or Roche Service representative.

COAGUCHEK, COBAS and COBAS H are trademarks of Roche.

All other trademarks are the property of their respective owners.

On the packaging, or the Handheld Base Unit you may encounter the following symbols, shown here with their meaning:

Caution, consult accompanying documents. Refer to safety-related notes in the instructions for use accompanying this product.

Temperature limitation (Store at)

Manufacturer

Date of manufacture

Catalog number

Global Trade Item Number

Unique device identifier

In vitro diagnostic medical device

Device for near-patient testing

Device not for self-testing

Complies with the provisions of the applicable EU regulations.

The Handheld Base Unit meets the requirements laid down in:

Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU.

Directive 2014/53/EU of the European Parliament and of the Council of 16 April 2014 on the harmonization of the laws of the Member States relating to the making available on the market of radio equipment and repealing Directive 1999/5/EC.

Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

Consult instructions for use

The system fulfills the Canadian and U.S. safety requirements in accordance with CAN/CSA C22.2 No. 61010-1 and UL 61010-1.

Eurasian Conformity. Demonstrates that the product meets the Eurasian Economic Union (EAEU)'s regulations and standards for customs clearance and trading.

12V  1.25A Power supply connection

This page intentionally left blank.

1	Introduction	9
	Purpose of document	9
	Intended Use	9
	Contents of the pack	9
	If you need help	9
	Labeling of safety information	10
	Safety information	11
	Disposal of the system	13
	General care	13
	Electrical safety	13
	Electromagnetic Interference	14
	HBU operating conditions	14
2	Overview of the Handheld Base Unit (HBU)	15
	Connection ports	17
	Configuration switch	18
3	Operating the device	19
3.1	Connecting the CoaguChek XS Plus, CoaguChek XS Pro, and CoaguChek Pro II meter	20
3.2	Connecting the cobas h 232 meter	21
	Important notes	23
	HBU self-check	24
4	Cleaning/Disinfection	25
	Introduction	25
	Cleaning/disinfecting agents	25
	Cleaning/disinfecting the Handheld Base Unit	26
5	Troubleshooting	27
6	Product specification	28
	Technical data and operating conditions	28
	Ordering	29
	Repairs	29
	Contact Roche	29
7	Information about software licenses	30

This page intentionally left blank.

1 Introduction

Purpose of document

This document describes the use of the Handheld Base Unit (also called **HBU** in this document).

Intended Use

The Handheld Base Unit is intended for storing, power supply, and data transfer to CoaguChek XS Plus/XS Pro/Pro II and **cobas h 232**. The Handheld Base Unit is intended for near-patient testing. Not for self-testing.

Contents of the pack

- Handheld Base Unit
- Operator's Manual
- USB data cable
- Extension piece for **cobas h 232** meter¹

If you need help

For all questions about the HBU that are not answered in this manual, contact your Roche representative. In order to expedite troubleshooting, please have your device, its serial number, and this manual ready when you call.

1. The extension piece for the **cobas h 232** meter is mandatory for meters with serial number ranges KS02XXXXX and KQ01XXXXX; optional for meters with serial number ranges KS04XXXXX and KQ03XXXXX

Labeling of safety information

This section explains how safety-related messages and information related to the proper handling of the system are presented in the Handheld Base Unit manual. Please read these passages carefully.

The safety alert symbol alone (without a signal word) promotes awareness to hazards which are generic or directs the reader to related safety information.

These symbols and signal words are used for specific hazards:

WARNING

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE

Indicates a hazardous situation which, if not avoided, may result in damage to the system.

Important information that is not safety relevant is presented against a colored background (without a symbol). Here you will find additional information on correct use of the Handheld Base Unit or useful tips.

Safety information

Operator qualification

Only trained healthcare professionals may operate the Handheld Base Unit. Operators must have received comprehensive instruction in the operation and care of the Handheld Base Unit.

Incident reporting

Inform your Roche representative and your local competent authority about any serious incidents which may occur when using this product.

WARNING

Protection against infection and blood-borne pathogens

Healthcare professionals using the Handheld Base Unit must be aware that any object coming into contact with human blood is a potential source of infection. Operators need to adhere to Standard Precautions when handling or using the Handheld Base Unit. All parts of this system should be considered potentially infectious and are capable of transmitting blood-borne pathogens between patients and between patients and healthcare professionals.

- Use gloves. Wear a new pair of clean gloves for each patient and control test.
 - Wash hands thoroughly with soap and water before putting on a new pair of gloves and performing the next patient test.
 - Follow all health and safety regulations in force locally.
-

 **CAUTION****Allergy or injury caused by reagents and other working solutions**

Direct contact with reagents, detergents, cleaning/disinfection solutions, or other working solutions may cause skin irritation or inflammation.

- Always use protective gloves.
- Observe the cautions given in the package inserts of the reagents and cleaning/disinfection solutions.
- If a reagent, control, linearity, or cleaning/disinfection solution comes into contact with your skin, wash it off immediately with water.
- Follow all health and safety regulations in force locally.

 **WARNING****Avoidance of electrical shock, fire, and explosions**

- Only use Roche Diagnostics original accessories (cables, power supply units, battery packs, and spare parts). Third-party cables, power supply units, and battery packs can cause the battery packs to explode or the devices to become damaged.
- Do not use loose power sockets or damaged power supply units, cables or plugs.
- Do not short circuit the power supply unit or the handheld base unit contacts.
- Do not drop the devices, the power supply units, or the battery packs and protect these against shaking and vibrations.

Disposal of the system

WARNING

Infection by a potentially biohazardous instrument

The Handheld Base Unit or its components must be treated as potentially biohazardous waste. Decontamination (i.e., a combination of processes including cleaning, disinfection and/or sterilization) is required before reuse, recycling, or disposal.

Dispose of the system or its components according to the appropriate local regulations. For more information, contact your Roche representative.

General care

NOTICE

Clean the Handheld Base Unit only with the solutions recommended (see chapter 4). Using other solutions may result in incorrect operation and possible system failure. Make sure that the Handheld Base Unit is thoroughly dried after cleaning or disinfecting.

Electrical safety

NOTICE

Never operate the Handheld Base Unit if the power supply unit or the attached cable are visibly damaged. If there is any visible damage, contact your local Roche service for inspections.

Electromagnetic Interference

The Handheld Base Unit fulfills the IEC 61326-2-6 requirements for emitted interference and interference immunity.

Do not use the Handheld Base Unit near strong electromagnetic fields, which could interfere with the proper operation of the Handheld Base Unit.

Electrostatic discharges may cause malfunction of the Handheld Base Unit.

HBU operating conditions

To ensure that the Handheld Base Unit functions properly, please observe the following guidelines:

- Only use the HBU at a room temperature between 12 °C and 32 °C (54 °F and 90 °F).
- Only use the device at a relative humidity between 10% and 85% (no condensation).
- When operating the device using the power supply unit, use only a voltage of 100V to 240 V ($\pm 10\%$), 50/60 Hz.

2 Overview of the Handheld Base Unit (HBU)

The HBU is used for:

- Charging the rechargeable meter battery pack when the meter is docked.
- Supporting communication with a data management system (DMS) via USB and Ethernet (RJ45) data ports.

The device has the following elements:

- 1 Charging contacts
- 2 Infrared window for communication with the meter
- 3 Status LED (lights up when power is connected)
- 4 Extension piece for **cobas h** 232 meter¹

HBU status LED

- Lights up red: Power supply is connected, application is starting up
- Lights up green: Ready
- Flashes red (once per second): Error
- Lights up blue: Configuration mode

1. The extension piece for the **cobas h** 232 meter is mandatory for meters with serial number ranges KS02XXXXX and KQ01XXXXX; optional for meters with serial number ranges KS04XXXXX and KQ03XXXXX

Connection ports

The following connection ports are located at the back of the Handheld Base Unit:

- 5** Power input jack for the power supply unit
- 6** Network connection - Ethernet/RJ45 port
- 7** USB port
- 8** Configuration switch

Configuration switch

The configuration switch on the back sets the HBU to one of three different modes of operation:


Pos	Mode
1	Configuration
2	Serial (ICI) & Ethernet (POCT1-A)
3	Serial (POCT1-A)

Always power off the HBU (remove the power cord) before changing the position of this switch.

For further information please refer to the *ROCHE HBU-BU-BUH Technical Note.pdf*, which can be found on the HBU (attached to a computer, acting as a removable disk). See page 19 for instructions on how to access the PDF file.

3 Operating the device

To help the System Administrator configure the HBU, technical documentation is provided. This technical note is stored on the HBU as a PDF file. You can access this file by connecting the HBU to a PC with the USB cable provided. Follow the directions below.


- 1 Make sure that the switch (8) is at position 1 (to your left in the illustration, see page 18). This is the default position when shipped.
- 2 Connect the power cord to the power input jack (5). After a maximum of 30 seconds the status LED changes to blue, indicating that the HBU is ready to be configured.
- 3 Connect the HBU to a PC using the USB cable. The HBU appears as a removable disk on the PC screen.
- 4 Copy the file *ROCHE HBU-BU-BUH Technical Note.pdf* to the desired location on your PC.
- 5 After changing the configuration settings, power off the HBU by removing the power cord from the input jack.
- 6 Move the configuration switch to either position 2 or 3, depending on your local settings. The new configuration will become active after powering on for the next time.
- 7 If you want the Handheld Base Unit to communicate
 - via USB: connect the USB cable connector to the appropriate port (7) or
 - via Ethernet: connect the RJ45 cable connector to the appropriate port (6).

You should plug in either a USB cable connector or an Ethernet cable connector, but not both at the same time.

Make sure the configuration switch is set to the correct position for the operating mode you selected.

3.1 Connecting the CoaguChek XS Plus, CoaguChek XS Pro, and CoaguChek Pro II meter

Place the CoaguChek XS Plus meter, the CoaguChek XS Pro meter or CoaguChek Pro II meter on the HBU, pushing it gently towards the charging contacts, until it snaps into position.


When the meter is correctly docked on the HBU and all cables are connected properly, battery charging and data transfer will start automatically.

A DMS connection has to be set up during initial HBU configuration in order to enable data transfer between meters and the DMS.

3.2 Connecting the cobas h 232 meter

For using the HBU with a **cobas h 232 meter**¹, first install the required extension piece.


- Install the HBU extension piece as shown above.

1. For meters with serial number ranges KS02XXXXX and KQ01XXXXX the extension piece has to be installed, otherwise the meter cannot be operated in the Handheld Base Unit.

For meters with serial number ranges KS04XXXXX and KQ03XXXXX the extension can, but does not have to be installed.

- Place the **cobas h 232** meter on the HBU, pushing it gently towards the charging contacts, until it snaps into position.

When the meter is correctly docked on the HBU and all cables are connected properly, battery charging and data transfer will start automatically.

A DMS connection has to be set up during initial HBU configuration in order to enable data transfer between meters and the DMS.

Important notes**Always...**

- ... handle the device carefully; avoid dropping or knocking it
- ... operate the device at an ambient temperature according to its operating conditions (see page 14 and/or page 28).
- ... place the device on a level, stable surface
- ... keep the device clean and dry

Never ...

- ... spill liquids on the device
- ... immerse the device in any liquid
- ... touch the charging contacts with metallic items, such as a paper clip, when plugged in
- ... move the HBU or undock the meter from the HBU while applying samples and performing a test.
- ... store the HBU at extreme temperatures.
- ... store the HBU in damp or humid conditions without protection.
- ... store the HBU in a dusty environment without protection.

HBU self-check

The Handheld Base Unit performs a self-check every 7 days (interval configurable, default value 7 days). The self-check causes the Handheld Base Unit software to reboot. To avoid interrupting communication between the meter and Handheld Base Unit, the software will only reboot when communication has been idle for at least 10 minutes. The test procedure itself takes about 3 minutes — in this time no communication between the meter and the host computer is possible. After completion of the self-check, the HBU will return automatically to the status it was in before the self-check began.

As part of the self-check, the integrity of the Handheld Base Unit's internal memory is checked. In case of a failure, the HBU will stay in USB mode even if the configuration switch is set to a different mode, and no serial or Ethernet communication of the Handheld Base Unit is possible.

4 Cleaning/Disinfection

Introduction

Healthcare professionals should wear gloves and follow their institution's infection control procedures when handling blood testing equipment. Use the following procedures as warranted.

NOTICE

Do not use any other cleaning/disinfecting agents (e.g. agents containing the chemicals ether, polyhexanide, or prepared solutions or wipes containing a mixture of bleach and detergent) on the Handheld Base Unit. Use of other cleaning/disinfecting agents could result in damage to the Handheld Base Unit.

Cleaning/disinfecting agents

Use only the following agents for cleaning/disinfecting the HBU:

Only for cleaning:

- Mild soapy water

For cleaning and disinfection:

- 70 % isopropyl alcohol
- Bacillo[®] AF
- Super Sani Cloth[®]

Use of other cleaning/disinfecting agents could result in damage to the HBU. Do not spray the HBU or into the connection ports.

Cleaning/disinfecting the Handheld Base Unit

- Unplug the HBU before cleaning/disinfecting.
- Wipe the surfaces with a soft cloth slightly dampened (**not wet**) with the cleaning/disinfecting agents as recommended above.
- Avoid wiping the electrical connections at the back of the base unit
- Ensure the HBU including all openings at the back (jack, ports) are **thoroughly dried** after cleaning or disinfecting. Visually verify that no solution is seen in the HBU power input jack, network connection port, USB port, and configuration switch at the completion of cleaning/disinfecting. If cleaning/disinfecting fluid is allowed to collect in any port, severe damage can occur to the meter and the HBU.
- If cleaning/disinfecting fluid does get into any opening (jack, ports), unplug the HBU, then **dry the components** with a dry cloth or gauze pad before returning the meter to the HBU.
- Plug in the HBU.

5 Troubleshooting

Display/symptom	Root cause	Possible solution
HBU light is not on	Not connected to power supply unit or power supply unit is defective, the HBU is defective, or mains power outlet not active.	Disconnect and connect the power supply again.
HBU light flashes	Charging contacts may be shorted.	Ensure that no object causes a possible short between the charging contacts. Disconnect and connect the power supply again.
Unit stays in USB configuration and LED flashes	The Handheld Base Unit performs a self-check regularly. In case of a failure, the HBU will stay in USB mode even if the configuration switch is set to a different mode, and no serial or Ethernet communication of the Handheld Base Unit is possible.	Check network connections, then disconnect and connect the power supply again.

If the error persists, contact your Roche representative.

6 Product specification

Technical data and operating conditions

Specification	Handheld Base Unit	Power supply unit
Height	82 mm / 3.22 in (max.)	35 mm / 1.38 in + AC plug (28-40 mm / 1.1-1.6 in) (max.)
Width	120 mm / 4.72 in (max.)	51 mm / 2.01 in (max.)
Length	228 mm / 8.97 in (max.) 295 mm / 11.61 in (max.) with extension	87 mm / 3.43 in (max.)
Weight	340 g / 11.99 oz (max.) 410 g / 14.46 oz (max.) with extension	170 g / 6 oz (max.)
User interface	LED	LED
Interfaces	Charging contacts IR port RJ45 Ethernet USB type B	DC connector Replaceable AC input contacts
Operating temperature	12 to 32 °C / 54 to 90 °F	0 to 40 °C / 32 to 104 °F
Storage temperature (long-term)	-5 to +45 °C / 23 to 113 °F	-5 to +45 °C / 23 to 113 °F
Humidity (operating)	10 to 85% (no condensation)	N/A
Relative humidity (storage)	< 95%	N/A
Input voltage	12 V DC	100 to 240 V AC
Input frequency	N/A	50 to 60 Hz
Input current	1.25 A (max.)	400-150 mA (max.) at 100-240 V AC / 50-60 Hz
Data transfer rate	IR: 9.6-115.2 Kbit/s Ethernet: 10/100 Mbit/s (auto-negotiate) full-duplex USB: 480 Mbit/s	N/A

Ordering

Below is a list of available accessories. For ordering information, contact your local sales representative (see section "Contact Roche").

Item	Description	REF
Handheld Base Unit		08408793XXX
Power Supply	Handheld base unit power supply Type: FW7555M/12	07006098001
	Handheld base unit power supply Type: FW8001M/12	08692432001

Repairs

Please note that repairs and other modifications to the meter may only be performed by persons authorized by Roche Diagnostics.

Contact Roche

For all questions about the Handheld Base Unit that are not answered in this manual, contact your Roche representative. A list of all Roche affiliates can be found at:
www.roche.com/about/business/roche_worldwide.htm

The Handheld Base Unit is manufactured for and distributed by:

In Australia:

Roche Diagnostics Australia Pty Limited
ABN 29 003 001 205
2 Julius Avenue
North Ryde, NSW, 2113

7 Information about software licenses

This product incorporates software modules developed under open source licenses. The source code of this software can be requested on a standard data exchange medium from the manufacturer at the following address:

Roche Diagnostics GmbH
Sandhofer Str. 116
68305 Mannheim
Germany

The complete license agreements are stored as a pdf file (file name "*license.pdf*") on the Handheld Base Unit. The file *license.pdf* is located in the same folder as the file *ROCHE HBU-BU-BUH Technical Note.pdf*.

You can access this file by connecting the Handheld Base Unit to a PC with the USB cable. For detailed instructions on how to do this, see the section "Operating the device" in this Operator's Manual.

Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

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